FLEX
Monoclonal Mouse
Anti-Thyroid
Transcription Factor, TTF-1
Clone 8G7G3/1
Ready-to-Use
(Dako Autostainer/Autostainer Plus)

Code IS056

English
Intended use

For in vitro diagnostic use.

FLEX Monoclonal Mouse Anti-Thyroid Transcription Factor, Clone 8G7G3/1 (1), Ready-to-Use, (Dako Autostainer/Autostainer Plus), is intended for use in immunohistochemistry together with Dako Autostainer/Autostainer Plus instruments. This antibody is useful for the identification of tumors from lung and thyroid. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary
and explanation

Thyroid transcription factor-1 (TTF-1) belongs to a family of homeodomain transcription factors and is selectively expressed in thyroid, lung and diencephalon. TTF-1 has been identified as a transcriptional regulator of thyroid-specific genes and has also been shown to be important in the activation of pulmonary-specific differentiation genes (2,3).

Refer to Dako's *General Instructions for Immunohistochemical Staining* or the detection system instructions of IHC procedures for: 1) Principle of Procedure, 2) Materials Required, Not Supplied, 3) Storage, 4) Specimen Preparation, 5) Staining Procedure, 6) Quality Control, 7) Troubleshooting, 8) Interpretation of Staining, 9) General Limitations.

Reagent provided

Ready-to-Use monoclonal mouse antibody provided in liquid form in a buffer containing stabilizing protein and 0.015 mol/L NaN₃.

Clone: 8G7G3/1 (1) Isotype: IgG1, kappa

Immunogen

Recombinant rat TTF-1 (1)

Specificity

Monoclonal anti-rat TTF-1 clone 8G7G3/1 (anti-TTF-1) has also been shown to react with rat, mouse and human TTF-1. In Western blotting, anti-TTF-1 specifically identified a 40 kDa band in immunoblots of nuclear extracts or whole cell lysates from the TTF-1 positive cell lines MLE 15 (mouse lung epithelial), H441-4 (human lung adenocarcinoma), H 345 (human small cell carcinoma) and rat type II pneumocyte cells. Cell lines lacking TTF-1 mRNA such as HeLa and 3T3 were unreactive. Anti-TTF-1 was also found to function in ELISA assays using recombinant TTF-1 (1).

Precautions

1. For professional users.
2. This product contains sodium azide (NaN₃), a chemical highly toxic in pure form. At product concentrations, though not classified as hazardous, sodium azide may react with lead and copper plumbing to form highly explosive build-ups of metal azides. Upon disposal, flush with large volumes of water to prevent metal azide build-up in plumbing.
3. As with any product derived from biological sources, proper handling procedures should be used.
4. Wear appropriate Personal Protective Equipment to avoid contact with eyes and skin.
5. Unused solution should be disposed of according to local, State and Federal regulations.

Storage

Store at 2-8 °C. Do not use after expiration date stamped on vial. If reagents are stored under any conditions other than those specified, the conditions must be verified by the user. There are no obvious signs to indicate instability of this product. Therefore, positive and negative controls should be run simultaneously with patient specimens. If unexpected staining is observed which cannot be explained by variations in laboratory procedures and a problem with the antibody is suspected, contact Dako Technical Support.

Specimen preparation
including materials
required but not
supplied

The antibody can be used for labeling formalin-fixed, paraffin-embedded tissue sections. Tissue specimens should be cut into sections of approximately 4 µm.

Pre-treatment with heat-induced epitope retrieval (HIER) is required. Optimal results are obtained by pretreating tissues using EnVision™ FLEX Target Retrieval Solution, High pH (10x) (Dako Autostainer/Autostainer Plus) (Code K8012/K8014).

Deparaffinized sections: Pre-treatment of deparaffinized formalin-fixed, paraffin-embedded tissue sections is recommended using Dako PT Link (Code PT100/PT101). For details, please refer to the PT Link User Guide.

Follow the pre-treatment procedure outlined in the package insert for EnVision™ FLEX Target Retrieval Solution, High pH (10x) (Dako Autostainer/Autostainer Plus) (Code K8012/K8014). The following parameters should be used for PT Link: Pre-heat temperature: 65 °C; epitope retrieval temperature and time: 97 °C for 20 (±1) minutes; cool down to 65 °C. Rinse sections with diluted room temperature EnVision™ FLEX Wash Buffer (10x) (Dako Autostainer/Autostainer Plus) (Code K8012).

Paraffin-embedded sections: As alternative specimen preparation, both deparaffinization and epitope retrieval can be performed in the PT Link using a modified procedure. See the PT Link User Guide for instructions. After the staining procedure has been completed, the sections must be dehydrated, cleared and mounted using permanent mounting medium.

The tissue sections should not dry out during the treatment or during the following immunohistochemical staining procedure. For greater adherence of tissue sections to glass slides, the use of Dako Silanized Slides (Code S3003) is recommended.

Staining procedure including materials required but not supplied

The recommended visualization system is EnVision™ FLEX+, Mouse, High pH (Dako Autostainer/Autostainer Plus) (Code K8012). The staining steps and incubation times are pre-programmed into the software of Dako Autostainer/Autostainer Plus instruments, using the following protocols:

Template protocol: FLEXRTU2 (200 uL dispense volume) or FLEXRTU3 (300 uL dispense volume)

Autoprogram (without counterstaining): TTF-1 or Autoprogram (with counterstaining): TTF-1H

The Auxiliary step should be set to "rinse buffer" in staining runs with ≤10 slides. For staining runs with >10 slides the Auxiliary step should be set to "none." This ascertains comparable wash times.

All incubation steps should be performed at room temperature. For details, please refer to the Operator's Manual for the dedicated instrument. If the protocols are not available on the used Dako Autostainer instrument, please contact Dako Technical Services.

Optimal conditions may vary depending on specimen and preparation methods, and should be determined by each individual laboratory. If the evaluating pathologist should desire a different staining intensity, a Dako Application Specialist/Technical Service Specialist can be contacted for information on re-programming of the protocol. Verify that the performance of the adjusted protocol is still valid by evaluating that the staining pattern is identical to the staining pattern described in "Performance characteristics."

Counterstaining in hematoxylin is recommended using EnVision™ FLEX Hematoxylin (Dako Autostainer/Autostainer Plus) (Code K8018). Non-aqueous, permanent mounting medium is recommended.

Positive and negative controls should be run simultaneously using the same protocol as the patient specimens. The positive control tissue should include thyroid follicular epithelial cells and the cells/structures should display reaction patterns as described for this tissue in "Performance characteristics" in all positive specimens. The recommended negative control reagent is FLEX Negative Control, Mouse (Dako Autostainer/Autostainer Plus) (Code IS750).

Staining interpretation

The cellular staining pattern is nuclear.

Performance characteristics

Normal tissues: Anti-TTF-1 has been shown to be immunoreactive with Type II cells and Clara cells of the lung and follicular cells from thyroid, but was unreactive with all other normal tissues examined. The following tissues were found to lack TTF-1 expression: pituitary, prostate, testes, adrenal gland, skin, breast, kidney, colon, liver, pancreas, small intestine, brain and stomach (1,4,5).

Abnormal tissues: TTF-1 expression has been demonstrated by IHC in neoplastic cells derived from tumors of the lung and thyroid. Anti-TTF-1 was found to react positively with the majority of pulmonary small cell carcinomas and primary and metastatic pulmonary adenocarcinomas, whereas a smaller proportion of large cell undifferentiated lung carcinomas (26%) were immunoreactive (1,5). TTF-1 has also been reported to react positively with squamous cell carcinomas of the lung (14%) (6). TTF-1 immunoreactivity was demonstrated in the majority of pulmonary atypical carcinoids but was rare in typical carcinoids (7). TTF-1 has also been detected in (3/3) thyroid papillary carcinomas (1). TTF-1 expression was not found in the majority of other tumors tested such as primary breast carcinomas, breast carcinomas metastatic to the lung, renal cell carcinomas, primary and metastatic colon and prostate adenocarcinomas, and malignant mesotheliomas (1,4,5,8,9). When using a polyclonal antibody to TTF-1, some focal immunoreactivity was observed in 1/66 gastric adenocarcinomas and 1/8 endometrial adenocarcinomas (8).

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References
Bibliographie
Literaturangaben

- Holzinger A, Dingle S, Bejarano PA, Miller M-A, Weaver TE, Di Lauro R, Whitsett JA. Monoclonal antibody to thyroid transcription factor-1: Production, characterization and usefulness in tumor diagnosis. Hybridoma 1996; 15(1):49
- Lazzaro D, Price M, De Felice M, Di Lauro R. The transcription factor TTF-1 is expressed at the onset of thyroid and lung morphogenesis and in restricted regions of the fetal brain. Development 1991; 113:1093
- Bohinski RJ, Di Lauro R, Whitsett JA. The lung-specific surfactant protein B gene promoter is a target for thyroid transcription factor-1 and hepatocyte nuclear factor-3, indicating common factors for organ-specific gene expression along the foregut axis. Mol Cell Biol 1994; 14(9):5671
- Bohinski RJ, Bejarano PA, Balko G, Warnick RE and Whitsett JA. Determination of lung as the primary site of cerebral metastatic adenocarcinomas using monoclonal antibody to thyroid transcription factor-1. J Neuro-Oncol. 1998; 40(3):227-31
- Khoor A, Whitsett JA, Stahlman MT, Stephenson M, Olson SJ, Cagle PT. Utility of surfactant protein B precursor and thyroid transcription factor 1 in differentiating adenocarcinoma of the lung from malignant mesothelioma. Hum Pathol 1999; 30(6):695-700
- Haque AK, Syed S, Lele SM, Freeman DH, Adegboyega PA. Immunohistochemical study of thyroid transcription factor-1 and HER2/neu in non-small cell lung cancer: strong thyroid transcription factor-1 expression predicts better survival. Appl. Immunohistochem Mol Morph 2002 June; 10(2):103
- Folpe AL, Hansen D, Gown AM, Schmidt RA. Expression of thyroid transcription factor-1 (TTF-1) is common in pulmonary atypical carcinoids and small cell carcinomas, but not typical carcinoids. Lab Invest 1998; 78(1):174A
- Bejarano PA, Baughman RP, Biddinger PW, Miller MA, Fenoglio-Preiser C, Al-Kafaji B, Di Lauro R, Whitsett JA. Surfactant proteins and thyroid transcription factor-1 in pulmonary and breast carcinomas. Mod Pathol 1996; 9(4):445
- Harlamert HA, Mira J, Yassin R, Bejarano PA, Baughman R, Miller MA, Whitsett J. Distinguishing primary lung adenocarcinomas from metastatic breast adenocarcinomas in cytology specimens using TTF-1 and cytokeratins 7 and 20. Mod Pathol 1997; 10(1):A185

 Catalogue number Référence du catalogue Bestellnummer	 Temperature limitation Limites de température Zulässiger Temperaturbereich	 In vitro diagnostic medical device Dispositif médical de diagnostic in vitro In-vitro-Diagnostikum	 Dako North America, Inc. 6392 Via Real Carpinteria, California 93013 USA Tel 805 566 6655 Fax 805 566 6688 Technical Support 800 424 0021 Customer Service 800 235 5763  Dako Denmark A/S Produktionsvej 42 DK-2600 Glostrup Denmark Tel +45 4485 9500 Fax +45 4485 9595 www.dako.com
 Manufacturer Fabricant Hersteller	 Batch code Code du lot Chargenbezeichnung	 Contains sufficient for <N> tests Contenu suffisant pour <n> tests Inhalt ausreichend für „n“ Ansätze	
 Use by Utiliser jusque Verwendbar bis	 Consult instructions for use Consulter les instructions d'utilisation Gebrauchsanweisung beachten	 Authorized representative in the European Community Représentant Autorisé dans la Communauté Européenne Autorisierter Repräsentant in der EU	

PT0020/ Rev C

Edition 09/07